

Table 1. Results of the Liquid Semen Dilutions

| Vasectomized Liquid Semen | | Dilution | Nonvasectomized Liquid Semen | |
|---------------------------|------------------------|------------------------------------------|------------------------------|------------------------|
| OneStep ABA Card® p30 | Seratec® PSA Semiquant | | OneStep ABA Card® p30 | Seratec® PSA Semiquant |
| + | + T > q | 1:100 | + | + T > q |
| + | + T > q | 1:1,000 | + | + T > q |
| + | + T > q | 1:10,000 | + | + T > q |
| + | + T > q | 1:50,000 | + | + T > q |
| + | + T > q | 1:100,000 | + | + T > q |
| + | + T > q | 1:200,000 | + | + T > q |
| + | + T > q | 1:400,000 | + | + T > q |
| + | + t < Q | 1:500,000 | - | + t < Q |
| - | + t < Q | 1:600,000 | - | + t < Q |
| - | + t < Q | 1:800,000 | - | + t < Q |
| - | + t < Q | 1:1,000,000 | - | + t < Q |
| - | - | Negative Control (sterile water only) | - | - |